

# FloorBed® 1





# FloorBed<sup>®</sup> 1

Reduce the risk of injury in the event of a patient fall

Reduce the risk of caregiver injury


Enable independence and mobility




## TECHNICAL SPECIFICATION

Dimensions:	925 x 2255 mm (36.4 x 88.8 in)	Lifting range:	71 - 650 mm (2.8 - 25.6 in) <i>(excl. mattress)</i>
Platform size:	900 x 2000 mm (35.4 x 78.7 in)	Profiling functions:	Backrest with Auto-regression, Kneebreak, Anti-Trendelenburg, Trendelenburg
Safe working load:	185 kg (29 st, 408 lbs)	Part number:	NSB-0-FL1-200
Patient weight:	150 kg (24 st, 330 lbs)		

## FloorBed® benefits:



Reduce bed fall to a roll with the aim to minimise injury



Enable independent leg lifting and reduce carer requirements



Promote independent bed transfers



## FloorBed® 1 PLUS

LIFTING RANGE	78 - 650 mm
MATTRESS SIZE	1050 x 2000 mm
MATTRESS PLATFORM EXTENSION	✓ + 200 mm
SIDE RAILS	✓
SAFE WORKING LOAD	260 kg (41 st)
PATIENT WEIGHT	225 kg (35 st)

## FloorBed® 1 Features



Width adjustable mattress guides



No trailing wires



Full nursing height



Concealed wheels



Backrest and kneebreak functions



Mattress platform height of just 7.1 cm – similar to a tennis ball!

## FloorBed® 1 Accessories



Integrated side rails made of soft material



Head and foot board bumpers may reduce impact and minimise the risk of injury



Safety mat, 22 cm high. Provides same height as a FloorBed® with 15 cm foam mattress



Safety sleeve may reduce risk of entrapment as the platform mechanism is protected



Bed lever to assist with bed mobility



Rotating bed lever for extra support in movement



Support bar to assist with bed mobility




Easy to transport


Visit our website at [www.floorbed.com](http://www.floorbed.com) for the full range of FloorBed® 1 accessories.

# Why FloorBed®?

## Nursing and Low



## FloorBed®



## When to use FloorBed®?

If the user:

- Is at risk of bed falls
- Has uncontrolled movements
- Has difficulties with leg lifting, difficulty with bed transfers and difficulty with bed mobility



## Injury risks when falling from bed

