



### **Description**

SOFT TOILET RAISER CONTACT PLUS NEO - 11 CM



### **Selling Points**

The Contact plus Néo toilet raisers brings a comfort never equalled on the market thanks to the Herdegen Softinject technology. With its double fixing lips, Contact Plus Néo quickly and securely fits on all standard toilet bowls to raise their height by 11 cm.

**UDI-DI : 3701127704426**

## WEIGHT - ECO-PERFORMANCE

Weight	2 kg
--------	------

## DIMENSIONS

Width	37 cm
Depth	41 cm
Height	16 cm
Dimensions (cm)	Toilet elevation : 11 cm

## MATERIAL & COATING

Seat	Injected HR foam
Material	PVC

## SAFETY CHARGE

Maximum patient weight	185 kg
------------------------	--------

## REGULATION

CE marking	CE
------------	----